



The Only Way to Fill



## DOUBLE PITCH FILLING AND STOPPERING MONOBLOCK MODEL « M53-4PP/2S »

**M.A.R.** calls attention of the visitors to the machine model **M53-4PP /2S/2A**, monobloc model among single or double pitch M53 machine series for dosage of **pharmaceutical substances in liquid form**.

Treatable vials formats from 2R to 100R (up to 250 ml vials upon appropriate modification).

The displayed version (double pitch) is equipped with 4 peristaltic pumps, but stainless steel or ceramic rotary piston volumetric dosing pumps or CIP SIP pumps can be installed on demand.

The machine is **equipped with 100% weight control, nitrogen stations, capping pressure control and a separation between production areas**.

The system is designed to accommodate an oRabs but has been designed also to be suitable for the installation under cRabs and Isolators.

Production speed 120 pieces per minute on 2R bottle.



**M.A.R. S.p.A.** | Viale Italia, 50 - 20020 Lainate (MI) Italy | Ph: + 39 029375371

E-Mail: [infosales@margroup.it](mailto:infosales@margroup.it) | [www.margroup.it](http://www.margroup.it)